



Installation instructions

System requirements

- IDEA3S supports Windows OS 2000 and XP.
- 50 MB free hard drive space

Installation

You may need to ask for administrator rights from your IT department before installing new programs on your PC

1. Double click on install_idea3s.zip or install_idea3s.exe on the download page and follow instructions
2. Open 'install_idea3s' folder on your hard drive
3. Click on setup.exe to start installation
4. You will be asked to agree or disagree to a disclaimer on the use of IDEA³S (if you disagree, you will not be able to continue with installing IDEA³S)
5. Follow instructions on screen to complete installation (accept all default names)

Start IDEA³S

1. Click on the IDEA³S icon on your desktop
2. You will be asked to provide a password; type "RMO" as the password

How to use IDEA³S

1. Go to the IDEA3S homepage and roll mouse over "About IDEA³S "
2. Choose "User's instruction" from drop down menu to view a PDF document titled "Help for using IDEA³S "

Notes for this version (v 1.0.8)

1. Homepage hyperlinks: Austin Health's antimicrobial-related guidelines are provided as examples under 'Hospital Guidelines'.
2. Internet connection is required to access:
 - Therapeutic Guidelines: Antibiotic 12th edition. At Austin Health, access to these guidelines are freely provided by the Clinicians Health Channel - health information website developed by the Dept. of Human Services, Victoria. Australia
 - Hyperlinks on the 'Creatinine clearance calculator' page – page is from Therapeutic Guidelines: Antibiotic 12th Edition
 - Hyperlinks on 'Intravenous Guidelines' page.
3. Patient numbers '123454', '123455' and '123456' are set up as test patients. To set up real patients, contact the developers.
4. Prescriber number '1234' is set up as a test doctor. To set up real doctors, contact the developers.

How to modify IDEA³S to suit your organisation

Refer to website for details on the IDEA³S training day. Alternatively, contact the developers for a copy of the IDEA³S Administrator Manual.

Updates

We will notify, by email, any updates to the application.

Contact Details

If you experience any problems or want to give feedback/comments, please don't hesitate to:

Email: idea3s.admin@austin.org.au or

Phone: 61.3.94965846